

Technical data

la belle 2

811/821


Compact hearing glasses module with 312 battery.
 For adjustment ranges, product details and available options, please consult our product information or visit our www.bruckhoff.com site

Notice: Measurements of digital hearing systems with pure tones can entail a wavy frequency response. This artefact is due to the use of a narrow band signal and does not reflect the actual achievement with wide-band input signals.


2 cm³ coupler data

ANSI/ASA S3.22.2014
 IEC 60118-7: 2015

Sound output level



Acoustic gain



		DIN EN 60118-0	DIN EN 60118-7
OSPL L _{Amax}	(db)	124	116
Gain max. (FOG 50)	(dB)	54	45
Equivalent input noise	(dB)	23	23
Frequency range	(Hz)	200-8000	200-8000
Battery current	(mA)	1.2	1.2